

ABX Minocal

LOT

CX 510

CAL



(Exp.) 2025-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2				9.4					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41				4.45					± 0.06
HB HGB	g/dl	13.1				13.2					± 0.2
	g/l	131				132					± 2
	mmol/l	8.14				8.20					± 0.12
HT HCT	%	37.0				37.4					± 1.0
	l/l	0.370				0.374					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242				249					± 10
VMP MPV	µm ³ ; fl	8.0				8.0					± 0.5

ABX Minocal

LOT

CX 510

CAL



(Exp.) 2025-10-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.42							± 0.06
HB HGB	g/dl		12.8							± 0.2
	g/l		128							± 2
	mmol/l		7.95							± 0.12
HT HCT	%		37.1							± 1.0
	l/l		0.371							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		250							± 10
VMP MPV *	µm ³ ; fl		8.0							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC