

ABX Minocal

LOT

CX 511

CAL



(Exp.) 2025-11-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8				8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.36				4.36					± 0.06
HB HGB	g/dl	13.0				12.9					± 0.2
	g/l	130				129					± 2
	mmol/l	8.07				8.01					± 0.12
HT HCT	%	36.6				36.4					± 1.0
	l/l	0.366				0.364					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220				241					± 10
VMP MPV	µm ³ ; fl	8.4				8.5					± 0.5

ABX Minocal

LOT

CX 511

CAL



(Exp.) 2025-11-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.34							± 0.06
HB HGB	g/dl		12.7							± 0.2
	g/l		127							± 2
	mmol/l		7.89							± 0.12
HT HCT	%		36.9							± 1.0
	l/l		0.369							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		234							± 10
VMP MPV *	µm ³ ; fl		8.0							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC