

ABX Minocal


LOT
CX 511

 Rev **1**
CAL

(Exp.) 2025-11-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0	9.0	8.9	9.1	8.9		8.9	8.9	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.37	4.32	4.39	4.41	4.36		4.43	4.39	± 0.06
HB HGB	g/dl	12.8	12.9	12.9	13.0	12.9	12.9		13.0	13.0	± 0.2
	g/l	128	129	129	130	129	129		130	130	± 2
	mmol/l	7.95	8.01	8.01	8.07	8.01	8.01		8.07	8.07	± 0.12
HT HCT	%	36.4	36.3	35.2	35.6	35.3	34.9		35.9	35.6	± 1.0
	l/l	0.364	0.363	0.352	0.356	0.353	0.349		0.359	0.356	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	235	228	231	244	230		232	231	± 10
VMP MPV	µm ³ ; fl	8.3	8.4	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.52 BACK / VERSO 1300090308