



Manufactured for
HORIBA Instruments Incorporated:
Pointe Brand
5449 Research Drive
Canton, MI 48188



LOT 522404

SAP 5390015524

REF T7535-CAL



2026-12-31

Calibrator Values*

	Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet
Calibrator A	0	mg/dL
	0.00	g/L
Calibrator B	79	mg/dL
	0.79	g/L
Calibrator C	158	mg/dL
	1.58	g/L
Calibrator D	264	mg/dL
	2.64	g/L
Calibrator E	370	mg/dL
	3.70	g/L
Calibrator F	528	mg/dL
	5.28	g/L

*Standardized against CAP/BCR/IFCC Reference Preparation for Proteins in Human Serum (RPPHS), lot 91/0619.