

Yumizen G FIB 5 ***for Yumizen G15X0x,*** ***Yumizen G8X0x and*** ***Yumizen G405***



Reference No.:		1300036384	
Lot No.:		950718 (rev. 1)	
Expiry date:		2027-07-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G15X0x/G8X0x/G405	g/L	sec
		4	9,7
		3	11,7
		2	15,3
		1	24,1
Dilution		10x	

Yumizen G FIB 5 ***for Yumizen G400*** ***DDi/G200***



Reference No.:		1300036384	
Lot No.:		950718 (rev. 1)	
Expiry date:		2027-07-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	9,6
		3	11,7
		2	15,5
		1	25,1
Dilution		10x	