

**LOT**

**CX 512**

**CAL**



**(Exp.) 2025-12-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500	H1500 H2500							
		HELO 1.x	HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.45</b>	<b>8.45</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.32</b>	<b>4.32</b>							± 0.06
HB HGB	g/dl	<b>12.7</b>	<b>12.7</b>							± 0.2
	g/l	<b>127</b>	<b>127</b>							± 2
	mmol/l	<b>7.89</b>	<b>7.89</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>230</b>	<b>230</b>							± 10
VGM MCV	µm <sup>3</sup> ; fl	<b>82.5</b>	<b>82.5</b>							± 4.0

Ref: TEMP-0387 Rev.53