

LOT

CX 512

CAL



(Exp.) 2025-12-05

(YYYY - MM - DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550		H500 OT	H500 CT H550	H560	H500 OT H500E OT	H500 CT H500E CT H550 H550E	H500 CRP	H560	
		v1.0 to v2.x			v3.x Only	v3.0 to v3.3		Since v4.x		Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75		8.60	8.60	8.60	8.65	8.65	8.65	8.65	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.38		4.33	4.33	4.33	4.34	4.34	4.34	4.34	± 0.06
HB HGB	g/dl	12.4		12.4	12.4	12.4	12.4	12.4	12.4	12.4	± 0.2
	g/l	124		124	124	124	124	124	124	124	± 2
	mmol/l	7.70		7.70	7.70	7.70	7.70	7.70	7.70	7.70	± 0.12
HT HCT	%	36.1		34.2	33.3	33.3	36.0	36.1	36.1	36.1	± 1.0
	l/l	0.361		0.342	0.333	0.333	0.360	0.361	0.361	0.361	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234		235	246	246	238	249	249	249	± 10
VMP MPV	μm ³ ; fl	9.3		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.53