

Yumizen G PT 5 for Yumizen G Line



Reference No.:		1300036338		
Lot No.:		950523 (rev. 1)		
Expiry date:		2027-05-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	13,1	
		50	18,6	
		33	24,3	
		25	29,6	
		1,19	MNPT 13,1	
Important notice				
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be improved by performing it at 37°C. 				