

Yumizen G DDi 2 ***for Yumizen G15X0x and*** ***Yumizen G405***



Reference No.:		1300036391	
Lot No.:		951003 (rev. 1)	
Expiry date:		2026-10-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G15X0x/G405	µg(FEU)/mL	ΔOD
		4.00	0.2775
		3.00	0.2060
		2.00	0.1270
		1.00	0.0575
		0.50	0.0280
		0.25	0.0145
Slope value		0,07	

Yumizen G DDi 2 ***for Yumizen G8X0x***



Reference No.:		1300036391	
Lot No.:		951003 (rev. 1)	
Expiry date:		2026-10-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G8X0x	µg(FEU)/mL	ΔOD
		4.00	0.2640
		3.00	0.2054
		2.00	0.1320
		1.00	0.0580
		0.50	0.0265
		0.25	0.0120
Slope value		0,06	



Yumizen G DDi 2 for Yumizen G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		951003 (rev. 1)	
Expiry date:		2026-10-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4.00	0.1225
		3.00	0.0925
		2.00	0.0565
		1.00	0.0255
		0.50	0.0132
		0.25	0.0076
		Slope value	0,02
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			