

**LOT**

**CX 513**

**CAL**



**(Exp.) 2026-01-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA	
		ES60 Care ST				60 8P 60 18P						
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.1</b>				<b>8.9</b>						± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.33</b>				<b>4.33</b>						± 0.06
HB HGB	g/dl	<b>13.1</b>				<b>13.0</b>						± 0.2
	g/l	<b>131</b>				<b>130</b>						± 2
	mmol/l	<b>8.14</b>				<b>8.07</b>						± 0.12
HT HCT	%	<b>36.8</b>				<b>36.8</b>						± 1.0
	l/l	<b>0.368</b>				<b>0.368</b>						± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>236</b>				<b>234</b>						± 10
VMP MPV	µm <sup>3</sup> ; fl	<b>8.2</b>				<b>7.5</b>						± 0.5

Ref: TEMP-0387 Rev.53