

LOT		CX 513		CAL		 (Exp.) 2026-01-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio							TOLERANCES TOLERANCE				
				MICROS		MICROS		PENTRA		PENTRA			PENTRA		
				EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60				60 60 C+ ES60	XL80	
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0	9.0	8.9	9.3	8.8		8.9	8.9	± 0.2			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.31	4.32	4.28	4.33	4.36	4.27		4.35	4.33	± 0.06			
HB	HGB	g/dl	12.8	12.8	12.8	13.0	12.9	12.8		13.0	13.0	± 0.2			
		g/l	128	128	128	130	129	128		130	130	± 2			
		mmol/l	7.95	7.95	7.95	8.07	8.01	7.95		8.07	8.07	± 0.12			
HT	HCT	%	36.2	36.1	34.5	35.5	35.1	34.4		35.7	35.5	± 1.0			
		l/l	0.362	0.361	0.345	0.355	0.351	0.344		0.357	0.355	± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	234	236	233	232	233	233		238	232	± 10			
VMP	MPV	µm ³ ; fl	8.4	8.3	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5			

Ref: TEMP-0387 Rev.53