

**LOT**

**CX 513**

Rev 1

**CAL**



**(Exp.) 2026-01-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 HELO 1.x	H1500 H2500 HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.35</b>	<b>8.35</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.36</b>	<b>4.36</b>							± 0.06
HB HGB	g/dl	<b>12.9</b>	<b>12.9</b>							± 0.2
	g/l	<b>129</b>	<b>129</b>							± 2
	mmol/l	<b>8.01</b>	<b>8.01</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>238</b>	<b>238</b>							± 10
VGM MCV	μm <sup>3</sup> ; fl	<b>82.5</b>	<b>82.5</b>							± 4.0

Ref: TEMP-0387 Rev.53