

LOT		CX 513		CAL		 (Exp.) 2026-01-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		YUMIZEN	YUMIZEN									
		H1500 H2500	H1500 H2500									
		HELO 1.x	HELO 2.0									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.35	8.35								± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.36	4.36								± 0.06	
HB HGB	g/dl	12.9	12.9								± 0.2	
	g/l	129	129								± 2	
	mmol/l	8.01	8.01								± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	238	238								± 10	
VGM MCV	µm ³ ; fl	82.5	82.5								± 4.0	

Ref. TEMP-0387 Rev.53