

LOT

CX 514

CAL



(Exp.) 2026-02-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.34				4.34					± 0.06
HB HGB	g/dl	12.8				12.8					± 0.2
	g/l	128				128					± 2
	mmol/l	7.95				7.95					± 0.12
HT HCT	%	36.0				36.0					± 1.0
	l/l	0.360				0.360					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237				235					± 10
VMP MPV	µm ³ ; fl	8.2				7.3					± 0.5

Ref: TEMP-0387 Rev.53