

## ABX Minocal

**LOT**

**CX 514**

**CAL**



**(Exp.) 2026-02-05**

Rev 1

(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.6</b>				<b>9.0</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.34</b>				<b>4.33</b>					± 0.06
HB HGB	g/dl	<b>12.7</b>				<b>12.7</b>					± 0.2
	g/l	<b>127</b>				<b>127</b>					± 2
	mmol/l	<b>7.89</b>				<b>7.89</b>					± 0.12
HT HCT	%	<b>35.6</b>				<b>35.9</b>					± 1.0
	l/l	<b>0.356</b>				<b>0.359</b>					± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>226</b>				<b>232</b>					± 10
VMP MPV	μm <sup>3</sup> ; fl	<b>8.2</b>				<b>8.4</b>					± 0.5

# ABX Minocal

**LOT**

**CX 514**

**CAL**



**(Exp.) 2026-02-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>8.6</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l		<b>4.35</b>							± 0.06
HB HGB	g/dl		<b>12.4</b>							± 0.2
	g/l		<b>124</b>							± 2
	mmol/l		<b>7.70</b>							± 0.12
HT HCT	%		<b>35.7</b>							± 1.0
	l/l		<b>0.357</b>							± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>228</b>							± 10
VMP MPV *	µm <sup>3</sup> ; fl		<b>7.9</b>							± 0.5

\* N/A for Micros Scil Vet ABC

Ref: TEMP-0387 Rev.53