

LOT

CX 514

Rev 1

CAL



(Exp.) 2026-02-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500	H1500 H2500							
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.35	8.35							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.36	4.36							± 0.06
HB HGB	g/dl	12.7	12.7							± 0.2
	g/l	127	127							± 2
	mmol/l	7.89	7.89							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	230							± 10
VGM MCV	µm ³ ; fl	82.0	82.0							± 4.0

Ref: TEMP-0387 Rev.53