

LOT		CX 514		CAL		(Exp.) 2026-02-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500 H2500	H1500 H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.35	8.35								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.36	4.36								± 0.06
HB HGB	g/dl	12.7	12.7								± 0.2
	g/l	127	127								± 2
	mmol/l	7.89	7.89								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	230								± 10
VGM MCV	μm ³ ; fl	82.0	82.0								± 4.0

Ref. TEMP-0387 Rev.53