


ABX Minotrol Retic

LOT RX 457

RETIC ADJUSTMENT TOOL

 (Exp.) **2026-03-05**
(YYYY - MM - DD)

Rev 1

PARAMETERS	UNITS	LEVEL 1				LEVEL 3			
		YUMIZEN	TOLERANCE		YUMIZEN	TOLERANCE			
		H2500 HELO 2.0			H2500 HELO 2.0				
RET	%	1.70			± 0.40			N/A	
MFI	-				N/A	0.335		± 0.030	