

LOT

CX 514

CAL



(Exp.) 2026-02-05

(YYYY - MM - DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550		H500 OT	H500 CT H550	H560	H500 OT H500E OT	H500 CT H500E CT H550 H550E	H500 CRP	H560	
		v1.0 to v2.x			v3.x Only	v3.0 to v3.3		Since v4.x		Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60		8.45	8.45	8.45	8.45	8.45	8.45	8.45	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42		4.30	4.30	4.30	4.35	4.35	4.35	4.35	± 0.06
HB HGB	g/dl	12.5		12.4	12.4	12.4	12.4	12.4	12.4	12.4	± 0.2
	g/l	125		124	124	124	124	124	124	124	± 2
	mmol/l	7.76		7.70	7.70	7.70	7.70	7.70	7.70	7.70	± 0.12
HT HCT	%	36.2		34.0	33.5	33.5	35.7	36.2	36.2	36.2	± 1.0
	l/l	0.362		0.340	0.335	0.335	0.357	0.362	0.362	0.362	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235		239	248	248	241	251	251	251	± 10
VMP MPV	μm ³ ; fl	9.9		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.53