



Yumizen G PT Liq 4 for Yumizen G Line



Reference No.:		1300036373		
Lot No.:		951014 (rev. 1)		
Expiry date:		2027-04-30		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	11,6	
		50	17,3	
		33	23,1	
		25	28,6	
		1,25	MNPT 11,6	
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				