

LOT

CX 515

Rev 1

CAL



(Exp.) 2026-03-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39				4.39					± 0.06
HB HGB	g/dl	12.6				12.6					± 0.2
	g/l	126				126					± 2
	mmol/l	7.82				7.82					± 0.12
HT HCT	%	36.2				36.2					± 1.0
	l/l	0.362				0.362					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	241				239					± 10
VMP MPV	µm ³ ; fl	8.2				7.3					± 0.5

Ref: TEMP-0387 Rev.53