

LOT		CX 515		CAL		 (Exp.) 2026-03-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio							TOLERANCES TOLERANCE		
				MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA			PENTRA
				EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60				60 60 C+ ES60
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1	9.1	9.1	9.1	9.0		9.0	9.1	± 0.2	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.37	4.34	4.33	4.34	4.42	4.35		4.40	4.34	± 0.06	
HB	HGB	g/dl	12.4	12.6	12.5	12.6	12.6	12.5		12.7	12.6	± 0.2	
		g/l	124	126	125	126	126	125		127	126	± 2	
		mmol/l	7.70	7.82	7.76	7.82	7.82	7.76		7.89	7.82	± 0.12	
HT	HCT	%	35.6	35.4	34.2	34.7	34.9	34.4		35.2	34.7	± 1.0	
		l/l	0.356	0.354	0.342	0.347	0.349	0.344		0.352	0.347	± 0.010	
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	239	245	236	231	244	237		239	231	± 10	
VMP	MPV	µm ³ ; fl	8.2	8.6	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

Ref: TEMP-0387 Rev.53