

LOT

CX 515

CAL



(Exp.) 2026-03-05

(YYYY - MM - DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550		H500 OT	H500 CT H550	H560	H500 OT H500E OT	H500 CT H500E CT H550 H550E	H500 CRP	H560	
		v1.0 to v2.x			v3.x Only	v3.0 to v3.3		Since v4.x		Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65		8.60	8.60	8.60	8.55	8.55	8.55	8.55	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39		4.35	4.35	4.35	4.42	4.37	4.37	4.37	± 0.06
HB HGB	g/dl	12.2		12.2	12.2	12.2	12.3	12.2	12.2	12.2	± 0.2
	g/l	122		122	122	122	123	122	122	122	± 2
	mmol/l	7.58		7.58	7.58	7.58	7.64	7.58	7.58	7.58	± 0.12
HT HCT	%	36.0		33.7	33.1	33.1	35.8	35.1	35.1	35.1	± 1.0
	l/l	0.360		0.337	0.331	0.331	0.358	0.351	0.351	0.351	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	249		254	257	257	251	257	257	257	± 10
VMP MPV	µm ³ ; fl	10.4		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.53