



**HORIBA ABX SAS**  
 Parc Euromédecine  
 Rue du Caducée  
 B.P. 7290  
 34184 Montpellier cedex 4  
 FRANCE  
 www.horiba-abx.com



A03A00012B

# Yumizen G APTT Liq 4 for YUMIZEN G Line



D 4 9 5 1 2 0 4 3 6 3 8 1 0 6 2 7 D



D 5 2 8 6 0 8 2 D

<b>Reference No.:</b>		<b>1300036381</b>
<b>Lot No.:</b>		<b>951204 (rev. 1)</b>
<b>Expiry date:</b>		<b>2027-06-30</b>
<b>MNAPTT</b>		
<b>Coagulometers</b>		<b>Sec</b>
<b>Optical</b>	<b>Yumizen G Line</b>	<b>28,6</b>
<b>Important notice</b>		
<p>The Mean Normal Activated Partial Thromboplastin Time (MNAPTT) is for information only. The value depends on the population, race, gender, sampling tube, etc. According to the CLSI every laboratory should determine its local MNAPTT.</p>		