

LOT		CX 516		CAL		 (Exp.) 2026-04-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500 H2500	H1500 H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.45	8.45								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42								± 0.06
HB HGB	g/dl	13.2	13.2								± 0.2
	g/l	132	132								± 2
	mmol/l	8.20	8.20								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	226	226								± 10
VGM MCV	µm ³ ; fl	83.0	83.0								± 4.0

Ref. TEMP-0387 Rev.53