

**LOT**

**CX 516**

Rev 1

**CAL**



**(Exp.) 2026-04-05**  
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500	H1500 H2500							
		HELO 1.x	HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.45</b>	<b>8.45</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.42</b>	<b>4.42</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.2</b>							± 0.2
	g/l	<b>132</b>	<b>132</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.20</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>226</b>	<b>226</b>							± 10
VGM MCV	µm <sup>3</sup> ; fl	<b>83.0</b>	<b>83.0</b>							± 4.0

Ref: TEMP-0387 Rev.53