

LOT

CX 517

CAL



(Exp.) 2026-05-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51				4.46					± 0.06
HB HGB	g/dl	13.2				13.2					± 0.2
	g/l	132				132					± 2
	mmol/l	8.20				8.20					± 0.12
HT HCT	%	37.4				37.2					± 1.0
	l/l	0.374				0.372					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	216				227					± 10
VMP MPV	μm ³ ; fl	8.3				8.3					± 0.5

ABX Minocal

LOT

CX 517

CAL



(Exp.) 2026-05-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.49							± 0.06
HB HGB	g/dl		13.1							± 0.2
	g/l		131							± 2
	mmol/l		8.14							± 0.12
HT HCT	%		37.3							± 1.0
	l/l		0.373							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		229							± 10
VMP MPV *	µm ³ ; fl		7.8							± 0.5

* N/A for Micros Scil Vet ABC

Ref: TEMP-0387 Rev.53