

LOT		CX 517		CAL		 (Exp.) 2026-05-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.45	8.45								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.53								± 0.06
HB HGB	g/dl	13.2	13.2								± 0.2
	g/l	132	132								± 2
	mmol/l	8.20	8.20								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	230								± 10
VGM MCV	µm ³ ; fl	82.0	82.0								± 4.0

Ref. TEMP-0387 Rev.53