



Manufactured for  
HORIBA Instruments Incorporated:  
Pointe Brand  
5449 Research Drive  
Canton, MI 48188



**LOT** 607601

**SAP** 5390015523



**REF** F7505-CTL



2027-11-30

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Level I Exp. Date: 2027-11-30 Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd			Level II Exp. Date: 2027-11-30 Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd		
		Mean	SD	Range	Mean	SD	Range
<b>Ferritin</b> (1st I.S. Standardization)	ng/mL	104	8.0	88-120	395	29.5	336-454
<b>Ferritin</b> (3rd I.S. Standardization)	ng/mL	124	9.5	105-143	470	35.3	400-541

Ferritin values are provided with two options: 1<sup>st</sup> or 3<sup>rd</sup> I.S. (International Standard) Values.

It is recommended to use the same standardized system as was chosen when selecting values for the Pointe Ferritin Calibrators (F7505-CAL).