

LOT

CX 518

CAL



(Exp.) 2026-06-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51				4.52					± 0.06
HB HGB	g/dl	13.5				13.5					± 0.2
	g/l	135				135					± 2
	mmol/l	8.38				8.38					± 0.12
HT HCT	%	37.9				38.4					± 1.0
	l/l	0.379				0.384					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	231				233					± 10
VMP MPV	µm ³ ; fl	8.2				7.5					± 0.5

Ref: TEMP-0387 Rev.53