

LOT

CX 518

CAL



(Exp.) 2026-06-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.2	9.1	9.4	9.1		9.1	9.1	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50	4.47	4.45	4.55	4.49		4.50	4.45	± 0.06
HB HGB	g/dl	13.3	13.2	13.4	13.4	13.5	13.3		13.5	13.4	± 0.2
	g/l	133	132	134	134	135	133		135	134	± 2
	mmol/l	8.26	8.20	8.32	8.32	8.38	8.26		8.38	8.32	± 0.12
HT HCT	%	37.8	37.8	36.7	36.3	36.9	36.4		36.7	36.3	± 1.0
	l/l	0.378	0.378	0.367	0.363	0.369	0.364		0.367	0.363	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	236	222	231	239	235		236	231	± 10
VMP MPV	µm ³ ; fl	8.4	8.4	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.53