

LOT		CX 518		CAL		 (Exp.) 2026-06-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500 H2500	H1500 H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65	8.65								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52								± 0.06
HB HGB	g/dl	13.4	13.4								± 0.2
	g/l	134	134								± 2
	mmol/l	8.32	8.32								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237	237								± 10
VGM MCV	µm ³ ; fl	82.5	82.5								± 4.0

Ref. TEMP-0387 Rev.53