



Yumizen G PT 5 for Yumizen G Line



D 4 9 5 1 2 0 1 3 6 3 3 8 1 2 2 7 D



D 5 1 0 0 1 2 6 5 0 0 1 8 4 3 3 0 D



D 6 2 4 4 2 5 0 3 0 0 1 2 6 1 1 9 D



D 7 1 4 9 D

Reference No.:		1300036338		
Lot No.:		951201 (rev. 1)		
Expiry date:		2027-12-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	12,6	
		50	18,4	
		33	24,4	
		25	30,0	
		1,19	MNPT 12,6	
Important notice				
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be improved by performing it at 37°C. 				