

Yumivet Cal

LOT VC 518 **CAL** (Exp.) 2026-06-05
 Rev 1 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Yumivet Cal		TOLERANCES TOLERANCE
		YUMIVET		
		VH2500		
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.65		± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52		± 0.06
HB HGB	g/dl	13.4		± 0.2
	g/l	134		± 2
	mmol/l	8.32		± 0.12
PLA PLT	10 ⁹ /mm ³ ;10 ⁹ /l	237		± 10
VGM MCV	μm ³ ; fl	82.5		± 4.0

Ref: TEMP-6057 Rev.01