



HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée
 BP 7290
 34184 Montpellier
 cedex 4
 FRANCE
 www.horiba-abx.com



A03A00010D

Yumizen G PT Reco 10 for Yumizen G Line



D 4 9 6 0 1 1 6 3 6 3 7 6 0 1 2 8 D



D 5 1 0 0 1 0 2 5 0 0 1 5 1 3 3 0 D



D 6 2 0 2 2 5 0 2 5 0 1 0 2 1 0 1 D



D 7 1 2 6 D

Reference No.:		1300036376		
Lot No.:		960116 (rev. 1)		
Expiry date:		2028-01-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	SEC
		100	10,2	
		50	15,1	
		33	20,2	
		25	25,0	
		1,01	MNPT	10,2
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p> <p>For Yumizen G200, Yumizen G400 and Yumizen G400 DDi, the Min. Time must be set to 7 seconds and the Lag time to 5 seconds in the PT setup. Please refer to the Yumizen G PT Reco 10 reagent notice.</p>				