

LOT

CX 519

CAL



(Exp.) 2026-07-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.31				4.36					± 0.06
HB HGB	g/dl	12.8				12.8					± 0.2
	g/l	128				128					± 2
	mmol/l	7.95				7.95					± 0.12
HT HCT	%	36.2				36.6					± 1.0
	l/l	0.362				0.366					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	231				232					± 10
VMP MPV	µm ³ ; fl	8.1				7.4					± 0.5

Ref: TEMP-0387 Rev.53