

LOT

CX 519

CAL



(Exp.) 2026-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500	H1500 H2500							
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.50	8.50							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.39							± 0.06
HB HGB	g/dl	12.7	12.7							± 0.2
	g/l	127	127							± 2
	mmol/l	7.89	7.89							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236	236							± 10
VGM MCV	µm ³ ; fl	81.5	81.5							± 4.0

Ref: TEMP-0387 Rev.53