



HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée
 B.P. 7290
 34184 Montpellier cedex 4
 FRANCE
 www.horiba-abx.com



A03A00009B

Yumizen G PT Liq 4 for Yumizen G Line



D 4 9 6 0 3 0 7 3 6 3 7 3 0 9 2 7 D



D 5 1 0 0 1 2 0 5 0 0 1 7 7 3 3 0 D



D 6 2 3 6 2 5 0 2 9 1 1 2 0 1 2 5 D



D 7 1 5 8 D

Reference No.:		1300036373		
Lot No.:		960307 (rev. 1)		
Expiry date:		2027-09-30		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,25	%
		100		12,0
		50		17,7
		33		23,6
		25		29,1
	MNPT	12,0		
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				