

LOT

CX 520

CAL



(Exp.) 2026-08-05

Rev **1**

(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52				4.52					± 0.06
HB HGB	g/dl	13.2				13.3					± 0.2
	g/l	132				133					± 2
	mmol/l	8.20				8.26					± 0.12
HT HCT	%	37.7				38.4					± 1.0
	l/l	0.377				0.384					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	223				225					± 10
VMP MPV	μm ³ ; fl	8.1				8.5					± 0.5

ABX Minocal

LOT

CX 520

CAL



(Exp.) 2026-08-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.5							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.53							± 0.06
HB HGB	g/dl		13.0							± 0.2
	g/l		130							± 2
	mmol/l		8.07							± 0.12
HT HCT	%		37.8							± 1.0
	l/l		0.378							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		226							± 10
VMP MPV *	µm ³ ; fl		7.7							± 0.5

Ref: TEMP-0387 Rev.53

*** N/A for Micros Scil Vet ABC**