


Yumivet Cal

LOT		VC 520	CAL		 (Exp.) 2026-08-05 <small>(YYYY - MM - DD)</small>
		Rev 1			
PARAMETRES PARAMETERS	UNITES UNITS	Yumivet Cal		TOLERANCES TOLERANCE	
		YUMIVET			
		VH2500			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8,30		± 0,20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4,54		± 0,06	
HB HGB	g/dl	13,1		± 0,2	
	g/l	131		± 2	
	mmol/l	8,14		± 0,12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	224		± 10	
VGM MCV	µm ³ ; fl	81,0		± 4,0	

Ref: TEMP-6057 Rev.01