


# ABX Minotrol Retic

**LOT** RX 459

**RETIC ADJUSTMENT TOOL**

 (Exp.) **2026-07-19**  
(YYYY - MM - DD)

Rev 2

PARAMETERS	UNITS	LEVEL 1				LEVEL 3			
		YUMIZEN			TOLERANCE	YUMIZEN			TOLERANCE
		H2500 HELO 2.0				H2500 HELO 2.0			
RET	%	<b>1.92</b>			± 0.40				N/A
MFI	-				N/A	<b>0.350</b>			± 0.030