

LAtrol™ Normal Control Plasma ***LAtrol®***

HORIBA Reference No.:	1300079435
Biomedica Reference No.:	816N
Lot No.:	260529
Expiry date:	2028/05/29

for Yumizen G800

Parameter	Unit	Range
DRVVT	sec	38.9 – 48.9
DRVVC	sec	34.8 – 42.0
Normalized Ratio	-	<1.16

for Yumizen G1500/G1550

Parameter	Unit	Range
DRVVT	sec	35.8 – 46.2
DRVVC	sec	33.1 – 44.9
Normalized Ratio	-	<1.16

Important notice

Each laboratory should establish its own DRVVT and DRVVC normal reference range (mean value $\pm$ 2 S.D.).

The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DRVVT (sample)} / \text{DRVVT (mean clot time test)}}{\text{DRVVC (sample)} / \text{DRVVC (mean clot time confirm)}}$$