

LOT		CX 521		CAL		(Exp.) 2026-09-05 (YYYY - MM - DD)					
Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58				4.57					± 0.06
HB HGB	g/dl	13.4				13.4					± 0.2
	g/l	134				134					± 2
	mmol/l	8.32				8.32					± 0.12
HT HCT	%	37.8				37.7					± 1.0
	l/l	0.378				0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	223				224					± 10
VMP MPV	µm ³ ; fl	7.8				7.0					± 0.5

Ref. TEMP-0387 Rev.53