

ABX Minocal

LOT		CX 521		CAL		(Exp.) 2026-09-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.5					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59				4.62					± 0.06
HB HGB	g/dl	13.4				13.3					± 0.2
	g/l	134				133					± 2
	mmol/l	8.32				8.26					± 0.12
HT HCT	%	37.9				37.9					± 1.0
	l/l	0.379				0.379					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	218				229					± 10
VMP MPV	µm ³ ; fl	7.7				8.0					± 0.5

ABX Minocal

LOT		CX 521		CAL		(Exp.) 2026-09-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.55								± 0.06
HB HGB	g/dl		13.2								± 0.2
	g/l		132								± 2
	mmol/l		8.20								± 0.12
HT HCT	%		38.0								± 1.0
	l/l		0.380								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		223								± 10
VMP MPV *	µm ³ ; fl		7.5								± 0.5

* N/A for Micros Scil Vet ABC

Ref: TEMP-0387 Rev.53