

| LOT                      |                                                       | CX 521      |         | CAL    |        | <div><div></div><div>(Exp.) 2026-09-05</div><div>( YYYY - MM - DD )</div></div> |        |        |                     |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|--------|---------------------|--------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                                                                                                                                                                    |        |        |                     |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA                                                                                                                                                             | PENTRA | PENTRA | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                                                                                                               | MS60   |        | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.2         | 9.1     | 9.1    | 9.0    | 9.5                                                                                                                                                                | 8.9    |        | 9.1                 | 9.0    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.55        | 4.56    | 4.47   | 4.56   | 4.58                                                                                                                                                               | 4.52   |        | 4.55                | 4.56   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.2        | 13.3    | 13.1   | 13.4   | 13.3                                                                                                                                                               | 13.2   |        | 13.4                | 13.4   | ± 0.2                   |
|                          | g/l                                                   | 132         | 133     | 131    | 134    | 133                                                                                                                                                                | 132    |        | 134                 | 134    | ± 2                     |
|                          | mmol/l                                                | 8.20        | 8.26    | 8.14   | 8.32   | 8.26                                                                                                                                                               | 8.20   |        | 8.32                | 8.32   | ± 0.12                  |
| HT HCT                   | %                                                     | 37.1        | 37.4    | 35.3   | 36.9   | 36.4                                                                                                                                                               | 35.9   |        | 36.6                | 36.9   | ± 1.0                   |
|                          | l/l                                                   | 0.371       | 0.374   | 0.353  | 0.369  | 0.364                                                                                                                                                              | 0.359  |        | 0.366               | 0.369  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 224         | 228     | 220    | 224    | 227                                                                                                                                                                | 226    |        | 226                 | 224    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 7.9         | 7.8     | N/A    | N/A    | N/A                                                                                                                                                                | N/A    |        | N/A                 | N/A    | ± 0.5                   |

Ref: TEMP-0387 Rev.53