

LOT		CX 521		CAL		 (Exp.) 2026-09-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500 H2500	H1500 H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.50	8.50								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57								± 0.06
HB HGB	g/dl	13.3	13.3								± 0.2
	g/l	133	133								± 2
	mmol/l	8.26	8.26								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	227								± 10
VGM MCV	µm ³ ; fl	82.0	82.0								± 4.0

Ref. TEMP-0387 Rev.53