



Yumizen G DDi 2 ***for Yumizen G15X0x and*** ***Yumizen G405***



D 4 9 6 0 3 3 9 3 6 3 9 1 0 3 2 7 D



D 5 4 0 0 2 1 8 0 3 0 0 1 8 1 0 2 D



D 6 0 0 1 3 0 0 1 0 0 6 5 0 0 5 D



D 7 0 0 3 3 8 0 2 5 0 1 9 5 1 7 3 D

Reference No.:		1300036391	
Lot No.:		960339 (rev. 1)	
Expiry date:		2027-03-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G15X0x/G405	4.00	0.2180
		3.00	0.1810
		2.00	0.1300
		1.00	0.0650
		0.50	0.0338
		0.25	0.0195
		Slope value	0,07

Yumizen G DDi 2 ***for Yumizen G8X0x***



D 4 9 6 0 3 3 9 3 6 3 9 1 0 3 2 7 D



D 5 4 0 0 2 2 2 3 0 0 1 7 5 0 2 D



D 6 0 0 1 2 2 0 1 0 0 6 4 0 0 5 D



D 7 0 0 3 4 2 0 2 5 0 1 9 0 1 6 3 D

Reference No.:		1300036391	
Lot No.:		960339 (rev. 1)	
Expiry date:		2027-03-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G8X0x	4.00	0.2222
		3.00	0.1750
		2.00	0.1220
		1.00	0.0640
		0.50	0.0342
		0.25	0.0190
		Slope value	0,06



Yumizen G DDi 2 ***for Yumizen*** ***G400 DDi/G200***

Reference No.:		1300036391	
Lot No.:		960339 (rev. 1)	
Expiry date:		2027-03-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4.00	0.0897
		3.00	0.0677
		2.00	0.0462
		1.00	0.0252
		0.50	0.0152
		0.25	0.0102
		Slope value	0,02
Important notice			
For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml . For more information, please contact your local HORIBA Medical representative.			