

## Section 1. Identification

**Product name** : Yumivet ReNew 0.5L  
**Product type** : Liquid.  
**Product code** : 1300159587  
**SAP Code** : 1300159587  
**Other means of identification** : Not available.  
**Product description** : 0.5 L

### Relevant identified uses of the substance or mixture and uses advised against

#### Identified uses

Cleaning solution intended for *in vitro* diagnostic use on HORIBA Veterinary blood cell counters.

#### Uses advised against

Not applicable.

**Supplier's details** : HORIBA ABX SAS  
 Parc Euromédecine  
 Rue du Caducée  
 BP 7290  
 34184 Montpellier Cedex 4  
 FRANCE  
 Tel: +33 (0) 4 67 14 15 16  
 Fax: +33 (0) 4 67 14 15 17

HORIBA Medical  
 9755 Research Drive  
 Irvine, CA 92618  
 Phone: + 1 888 903 5001  
 Fax: +1 949 753 0533

**e-mail address of person responsible for this SDS** : documentation.med@horiba.com  
 hotline-med.us@horiba.com

### National advisory body/Poison Center

**Emergency telephone number** : +1 800 424 9300 (CHEMTREC®)

#### Supplier

**Telephone number** : + 800 67 14 15 16

## Section 2. Hazards identification

**OSHA/HCS status** : This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).

**Classification of the substance or mixture** : CORROSIVE TO METALS - Category 1  
 AQUATIC HAZARD (LONG-TERM) - Category 3

### GHS label elements

## Section 2. Hazards identification

**Hazard pictograms** :



**Signal word** : Warning

**Hazard statements** : May be corrosive to metals.  
Harmful to aquatic life with long lasting effects.

### Precautionary statements

**Prevention** : Keep only in original packaging. Avoid release to the environment.

**Response** : Absorb spillage to prevent material damage.

**Storage** : Store in a corrosion resistant container with a resistant inner liner.

**Disposal** : Not applicable.

**Hazards not otherwise classified** : None known.

## Section 3. Composition/information on ingredients

**Substance/mixture** : Mixture

**Other means of identification** : Not available.

| Ingredient name                        | %  | CAS number |
|----------------------------------------|----|------------|
| sodium hypochlorite solution Cl active | <5 | 7681-52-9  |
| sodium hydroxide                       | ≤1 | 1310-73-2  |

Any concentration shown as a range is to protect confidentiality or is due to batch variation.

There are no additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified and hence require reporting in this section.

Occupational exposure limits, if available, are listed in Section 8.

## Section 4. First aid measures

### Description of necessary first aid measures

**Eye contact** : Immediately flush eyes with plenty of water, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical attention if irritation occurs.

**Inhalation** : Remove victim to fresh air and keep at rest in a position comfortable for breathing.

**Skin contact** : Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Get medical attention if symptoms occur.

**Ingestion** : Wash out mouth with water. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Do not induce vomiting unless directed to do so by medical personnel.

### Most important symptoms/effects, acute and delayed

#### Potential acute health effects

**Eye contact** : No known significant effects or critical hazards.

**Inhalation** : No known significant effects or critical hazards.

**Skin contact** : No known significant effects or critical hazards.

**Ingestion** : No known significant effects or critical hazards.

#### Over-exposure signs/symptoms

## Section 4. First aid measures

|                     |                     |
|---------------------|---------------------|
| <b>Eye contact</b>  | : No specific data. |
| <b>Inhalation</b>   | : No specific data. |
| <b>Skin contact</b> | : No specific data. |
| <b>Ingestion</b>    | : No specific data. |

### Indication of immediate medical attention and special treatment needed, if necessary

|                                   |                                                                                                                             |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Notes to physician</b>         | : Treat symptomatically. Contact poison treatment specialist immediately if large quantities have been ingested or inhaled. |
| <b>Specific treatments</b>        | : No specific treatment.                                                                                                    |
| <b>Protection of first-aiders</b> | : No action shall be taken involving any personal risk or without suitable training.                                        |

See toxicological information (Section 11)

## Section 5. Fire-fighting measures

### Extinguishing media

|                                       |                                                                 |
|---------------------------------------|-----------------------------------------------------------------|
| <b>Suitable extinguishing media</b>   | : Use an extinguishing agent suitable for the surrounding fire. |
| <b>Unsuitable extinguishing media</b> | : None known.                                                   |

|                                                   |                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Specific hazards arising from the chemical</b> | : In a fire or if heated, a pressure increase will occur and the container may burst. This material is harmful to aquatic life with long lasting effects. Fire water contaminated with this material must be contained and prevented from being discharged to any waterway, sewer or drain. |
|---------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                                                 |                                                                                                              |
|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Hazardous thermal decomposition products</b> | : Decomposition products may include the following materials:<br>halogenated compounds<br>metal oxide/oxides |
|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------|

|                                                     |                                                                                                                                                                                               |
|-----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Special protective actions for fire-fighters</b> | : Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable training. |
|-----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                                                       |                                                                                                                                                                       |
|-------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Special protective equipment for fire-fighters</b> | : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode. |
|-------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                  |
|---------------|------------------|
| <b>Remark</b> | : Non-flammable. |
|---------------|------------------|

## Section 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

|                                    |                                                                                                                                                                                                                                                                             |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>For non-emergency personnel</b> | : No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Put on appropriate personal protective equipment. |
| <b>For emergency responders</b>    | : If specialized clothing is required to deal with the spillage, take note of any information in Section 8 on suitable and unsuitable materials. See also the information in "For non-emergency personnel".                                                                 |

|                                  |                                                                                                                                                                                                                                                                                                                        |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Environmental precautions</b> | : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air). Water polluting material. May be harmful to the environment if released in large quantities. |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### Methods and materials for containment and cleaning up

## Section 6. Accidental release measures

- Small spill** : Stop leak if without risk. Move containers from spill area. Dilute with water and mop up if water-soluble. Alternatively, or if water-insoluble, absorb with an inert dry material and place in an appropriate waste disposal container. Absorb spillage to prevent material damage. Dispose of via a licensed waste disposal contractor.
- Large spill** : Stop leak if without risk. Move containers from spill area. Absorb spillage to prevent material damage. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows. Contain and collect spillage with non-combustible, absorbent material e.g. sand, earth, vermiculite or diatomaceous earth and place in container for disposal according to local regulations (see Section 13). Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. Note: see Section 1 for emergency contact information and Section 13 for waste disposal.

## Section 7. Handling and storage

### Precautions for safe handling

- Protective measures** : Put on appropriate personal protective equipment (see Section 8). Do not ingest. Avoid contact with eyes, skin and clothing. Avoid breathing vapor or mist. Avoid release to the environment. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Keep away from acids. Empty containers retain product residue and can be hazardous. Do not reuse container. Absorb spillage to prevent material damage.

- Advice on general occupational hygiene** : Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Remove contaminated clothing and protective equipment before entering eating areas. See also Section 8 for additional information on hygiene measures.

- Conditions for safe storage, including any incompatibilities** : Store between the following temperatures: 18 to 25°C (64,4 to 77°F). Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see Section 10) and food and drink. Store in a corrosion resistant container with a resistant inner liner. Separate from acids. Keep away from metals. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination. See Section 10 for incompatible materials before handling or use.

## Section 8. Exposure controls/personal protection

### Control parameters

#### Occupational exposure limits

| Ingredient name                                                | Exposure limits                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| sodium hypochlorite solution Cl active<br><br>sodium hydroxide | <b>OARS WEEL (United States, 4/2022).</b><br>STEL: 2 mg/m <sup>3</sup> 15 minutes.<br><b>ACGIH TLV (United States, 1/2023).</b><br>C: 2 mg/m <sup>3</sup><br><b>OSHA PEL 1989 (United States, 3/1989).</b><br>CEIL: 2 mg/m <sup>3</sup><br><b>NIOSH REL (United States, 10/2020).</b><br>CEIL: 2 mg/m <sup>3</sup><br><b>OSHA PEL (United States, 5/2018).</b><br>TWA: 2 mg/m <sup>3</sup> 8 hours.<br><b>CAL OSHA PEL (United States, 5/2018).</b><br>C: 2 mg/m <sup>3</sup> |

### Biological exposure indices

## Section 8. Exposure controls/personal protection

No exposure indices known.

|                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Appropriate engineering controls</b>      | : Good general ventilation should be sufficient to control worker exposure to airborne contaminants.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Environmental exposure controls</b>       | : Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.                                                                                                                                                                                                                                                                                                                                                    |
| <b><u>Individual protection measures</u></b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Hygiene measures</b>                      | : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.                                                                                                                                                                                                                                                                                     |
| <b>Eye/face protection</b>                   | : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts. If contact is possible, the following protection should be worn, unless the assessment indicates a higher degree of protection: safety glasses with side-shields. Recommended: Tightly sealed goggles according to EN 166                                                                                                                                                                                                                                                             |
| <b><u>Skin protection</u></b>                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Hand protection</b>                       | : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary. Considering the parameters specified by the glove manufacturer, check during use that the gloves are still retaining their protective properties. It should be noted that the time to breakthrough for any glove material may be different for different glove manufacturers. In the case of mixtures, consisting of several substances, the protection time of the gloves cannot be accurately estimated. > 8 hours (breakthrough time): Protective gloves according to EN 374 |
| <b>Body protection</b>                       | : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Other skin protection</b>                 | : Appropriate footwear and any additional skin protection measures should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Respiratory protection</b>                | : Based on the hazard and potential for exposure, select a respirator that meets the appropriate standard or certification. Respirators must be used according to a respiratory protection program to ensure proper fitting, training, and other important aspects of use. Recommended: Combination filtering device according to EN 14387                                                                                                                                                                                                                                                                                                                                |

## Section 9. Physical and chemical properties and safety characteristics

The conditions of measurement of all properties are at standard temperature and pressure unless otherwise indicated.

### Appearance

|                                     |                  |
|-------------------------------------|------------------|
| <b>Physical state</b>               | : Liquid.        |
| <b>Color</b>                        | : Yellowish.     |
| <b>Odor</b>                         | : Chlorine       |
| <b>Odor threshold</b>               | : Not available. |
| <b>pH</b>                           | : 12 to 13       |
| <b>Melting point/freezing point</b> | : Not available. |

## Section 9. Physical and chemical properties and safety characteristics

|                                                                |                                                                         |
|----------------------------------------------------------------|-------------------------------------------------------------------------|
| <b>Boiling point, initial boiling point, and boiling range</b> | : Not available.                                                        |
| <b>Flash point</b>                                             | : Not applicable.                                                       |
| <b>Flammability</b>                                            | : Non-flammable.                                                        |
| <b>Lower and upper explosion limit/flammability limit</b>      | : Not available.                                                        |
| <b>Vapor pressure</b>                                          | : Not available.                                                        |
| <b>Relative vapor density</b>                                  | : Not available.                                                        |
| <b>Relative density</b>                                        | : 1,02                                                                  |
| <b>Density</b>                                                 | : 1,02 g/cm <sup>3</sup> [20°C (68°F)]                                  |
| <b>Solubility in water</b>                                     | : Not available.                                                        |
| <b>Miscible with water</b>                                     | : Yes.                                                                  |
| <b>Partition coefficient: n-octanol/water</b>                  | : Not applicable.                                                       |
| <b>Auto-ignition temperature</b>                               | : Not applicable.                                                       |
| <b>Decomposition temperature</b>                               | : Not available.                                                        |
| <b>Viscosity</b>                                               | : Not available.                                                        |
| <b>Corrosion</b>                                               | : 56,7 [aluminum] mm per year (55°C)<br>5,3 [Steel.] mm per year (55°C) |

### Particle characteristics

|                             |                   |
|-----------------------------|-------------------|
| <b>Median particle size</b> | : Not applicable. |
|-----------------------------|-------------------|

## Section 10. Stability and reactivity

|                                           |                                                                                                        |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Reactivity</b>                         | : No specific test data related to reactivity available for this product or its ingredients.           |
| <b>Chemical stability</b>                 | : The product is stable.                                                                               |
| <b>Possibility of hazardous reactions</b> | : Under normal conditions of storage and use, hazardous reactions will not occur.                      |
| <b>Conditions to avoid</b>                | : No specific data.                                                                                    |
| <b>Incompatible materials</b>             | : Reactive or incompatible with the following materials:<br>acids<br>metals                            |
| <b>Other information</b>                  | : Yumivet BasoCare                                                                                     |
| <b>Hazardous decomposition products</b>   | : Under normal conditions of storage and use, hazardous decomposition products should not be produced. |

## Section 11. Toxicological information

### Information on toxicological effects

#### Acute toxicity

Not available.

#### Irritation/Corrosion

## Section 11. Toxicological information

| Product/ingredient name                                                                     | Result                             | Species | Score | Exposure         | Observation |
|---------------------------------------------------------------------------------------------|------------------------------------|---------|-------|------------------|-------------|
| Yumivet ReNew 0.5L<br><br>sodium hypochlorite solution<br>Cl active<br><br>sodium hydroxide | Eyes - Cornea opacity              | Rabbit  | 0,23  | -                | -           |
|                                                                                             | Eyes - Edema of the conjunctivae   | Rabbit  | 0,53  | -                | -           |
|                                                                                             | Eyes - Iris lesion                 | Rabbit  | 0     | -                | -           |
|                                                                                             | Eyes - Redness of the conjunctivae | Rabbit  | 1,53  | -                | -           |
|                                                                                             | Skin - Edema                       | Rabbit  | 0     | -                | -           |
|                                                                                             | Skin - Erythema/Eschar             | Rabbit  | 1,8   | -                | -           |
|                                                                                             | Eyes - Mild irritant               | Rabbit  | -     | 1.31 mg          | -           |
|                                                                                             | Eyes - Moderate irritant           | Rabbit  | -     | 10 mg            | -           |
|                                                                                             | Eyes - Mild irritant               | Rabbit  | -     | 400 ug           | -           |
|                                                                                             | Eyes - Severe irritant             | Monkey  | -     | 24 hours 1 %     | -           |
|                                                                                             | Eyes - Severe irritant             | Rabbit  | -     | 1 %              | -           |
|                                                                                             | Eyes - Severe irritant             | Rabbit  | -     | 0,5 minutes 1 mg | -           |
|                                                                                             | Eyes - Severe irritant             | Rabbit  | -     | 24 hours 50 ug   | -           |
|                                                                                             | Skin - Mild irritant               | Human   | -     | 24 hours 2 %     | -           |
|                                                                                             | Skin - Severe irritant             | Rabbit  | -     | 24 hours 500 mg  | -           |

### Sensitization

Not available.

### Mutagenicity

Not available.

### Carcinogenicity

Not available.

### Classification

| Product/ingredient name                   | OSHA | IARC | NTP |
|-------------------------------------------|------|------|-----|
| sodium hypochlorite solution<br>Cl active | -    | 3    | -   |

### Reproductive toxicity

Not available.

### Teratogenicity

Not available.

### Specific target organ toxicity (single exposure)

| Product/ingredient name                | Category   | Route of exposure | Target organs                |
|----------------------------------------|------------|-------------------|------------------------------|
| sodium hypochlorite solution Cl active | Category 3 | -                 | Respiratory tract irritation |

### Specific target organ toxicity (repeated exposure)

Not available.

### Aspiration hazard

Not available.

Information on the likely routes of exposure : Not available.

## Section 11. Toxicological information

### Potential acute health effects

|                     |                                                     |
|---------------------|-----------------------------------------------------|
| <b>Eye contact</b>  | : No known significant effects or critical hazards. |
| <b>Inhalation</b>   | : No known significant effects or critical hazards. |
| <b>Skin contact</b> | : No known significant effects or critical hazards. |
| <b>Ingestion</b>    | : No known significant effects or critical hazards. |

### Symptoms related to the physical, chemical and toxicological characteristics

|                     |                     |
|---------------------|---------------------|
| <b>Eye contact</b>  | : No specific data. |
| <b>Inhalation</b>   | : No specific data. |
| <b>Skin contact</b> | : No specific data. |
| <b>Ingestion</b>    | : No specific data. |

### Delayed and immediate effects and also chronic effects from short and long term exposure

#### Short term exposure

|                                    |                  |
|------------------------------------|------------------|
| <b>Potential immediate effects</b> | : Not available. |
| <b>Potential delayed effects</b>   | : Not available. |

#### Long term exposure

|                                    |                  |
|------------------------------------|------------------|
| <b>Potential immediate effects</b> | : Not available. |
| <b>Potential delayed effects</b>   | : Not available. |

### Potential chronic health effects

Not available.

|                              |                                                     |
|------------------------------|-----------------------------------------------------|
| <b>General</b>               | : No known significant effects or critical hazards. |
| <b>Carcinogenicity</b>       | : No known significant effects or critical hazards. |
| <b>Mutagenicity</b>          | : No known significant effects or critical hazards. |
| <b>Reproductive toxicity</b> | : No known significant effects or critical hazards. |

### Numerical measures of toxicity

#### Acute toxicity estimates

N/A

## Section 12. Ecological information

### Toxicity

| Product/ingredient name                   | Result                             | Species                                                                       | Exposure |
|-------------------------------------------|------------------------------------|-------------------------------------------------------------------------------|----------|
| sodium hypochlorite solution<br>Cl active | Acute EC50 0,67 mg/l Marine water  | Algae - Phaeodactylum<br>tricornutum - Exponential growth<br>phase            | 96 hours |
|                                           | Acute EC50 0,01 mg/l Fresh water   | Daphnia - Daphnia magna -<br>Embryo                                           | 48 hours |
|                                           | Acute LC50 56,4 mg/l Marine water  | Crustaceans - Palaemonetes<br>pugio                                           | 48 hours |
|                                           | Acute LC50 32 µg/l Marine water    | Fish - Oncorhynchus kisutch -<br>Juvenile (Fledgling, Hatchling,<br>Weanling) | 96 hours |
|                                           | Chronic NOEC 0,5 mg/l Marine water | Algae - Isochrysis galbana -<br>Exponential growth phase                      | 96 hours |

Date of issue/Date of revision : 10/03/2026

8/12

## Section 12. Ecological information

|                  |                                                                    |                                                                   |                     |
|------------------|--------------------------------------------------------------------|-------------------------------------------------------------------|---------------------|
| sodium hydroxide | Chronic NOEC 0,1 ppm Fresh water<br>Acute LC50 125 ppm Fresh water | Fish - Cyprinus carpio - Young<br>Fish - Gambusia affinis - Adult | 30 days<br>96 hours |
|------------------|--------------------------------------------------------------------|-------------------------------------------------------------------|---------------------|

### Persistence and degradability

Not available.

### Bioaccumulative potential

Not available.

### Mobility in soil

**Soil/water partition coefficient (K<sub>oc</sub>)** : Not available.

**Other adverse effects** : No known significant effects or critical hazards.

## Section 13. Disposal considerations

**Disposal methods** : The generation of waste should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Waste should not be disposed of untreated to the sewer unless fully compliant with the requirements of all authorities with jurisdiction. Waste packaging should be recycled. Incineration or landfill should only be considered when recycling is not feasible. This material and its container must be disposed of in a safe way. Care should be taken when handling emptied containers that have not been cleaned or rinsed out. Empty containers or liners may retain some product residues. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

## Section 14. Transport information

|                                   | <b>DOT<br/>Classification</b>                                                                | <b>TDG<br/>Classification</b>                                                                | <b>Mexico<br/>Classification</b>                                                             | <b>IMDG</b>                                                                                  | <b>IATA</b>                                                                                  |
|-----------------------------------|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
| <b>UN number</b>                  | UN3266                                                                                       | UN3266                                                                                       | UN3266                                                                                       | UN3266                                                                                       | UN3266                                                                                       |
| <b>UN proper shipping name</b>    | CORROSIVE LIQUID, BASIC, INORGANIC, N.O. S (sodium hypochlorite, solution, sodium hydroxide) | CORROSIVE LIQUID, BASIC, INORGANIC, N.O. S (sodium hypochlorite, solution, sodium hydroxide) | CORROSIVE LIQUID, BASIC, INORGANIC, N.O. S (sodium hypochlorite, solution, sodium hydroxide) | CORROSIVE LIQUID, BASIC, INORGANIC, N.O. S (sodium hypochlorite, solution, sodium hydroxide) | CORROSIVE LIQUID, BASIC, INORGANIC, N.O. S (sodium hypochlorite, solution, sodium hydroxide) |
| <b>Transport hazard class(es)</b> | 8<br>     | 8<br>     | 8<br>     | 8<br>   | 8<br>   |
| <b>Packing group</b>              | III                                                                                          | III                                                                                          | III                                                                                          | III                                                                                          | III                                                                                          |
| <b>Environmental hazards</b>      | No.                                                                                          | No.                                                                                          | No.                                                                                          | No.                                                                                          | No.                                                                                          |

### Additional information

## Section 14. Transport information

- DOT Classification** : **Reportable quantity** 5984,8 lbs / 2717,1 kg [703,71 gal / 2663,8 L]. Package sizes shipped in quantities less than the product reportable quantity are not subject to the RQ (reportable quantity) transportation requirements.
- TDG Classification** : Product classified as per the following sections of the Transportation of Dangerous Goods Regulations: 2.40-2.42 (Class 8).
- IMDG** : **Remarks** Limited quantity

**Special precautions for user** : **Transport within user's premises:** always transport in closed containers that are upright and secure. Ensure that persons transporting the product know what to do in the event of an accident or spillage.

**Transport in bulk according to IMO instruments** : Not available.

## Section 15. Regulatory information

**U.S. Federal regulations** : **TSCA 8(a) CDR Exempt/Partial exemption:** Not determined  
**Clean Water Act (CWA) 311:** sodium hypochlorite, solution; sodium hydroxide

**Clean Air Act Section 112 (b) Hazardous Air Pollutants (HAPs)** : Not listed

**Clean Air Act Section 602 Class I Substances** : Not listed

**Clean Air Act Section 602 Class II Substances** : Not listed

**DEA List I Chemicals (Precursor Chemicals)** : Not listed

**DEA List II Chemicals (Essential Chemicals)** : Not listed

### SARA 302/304

#### Composition/information on ingredients

No products were found.

**SARA 304 RQ** : Not applicable.

### SARA 311/312

**Classification** : CORROSIVE TO METALS - Category 1

#### Composition/information on ingredients

| Name                                   | %  | Classification                                                                                                                                                                                      |
|----------------------------------------|----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| sodium hypochlorite solution Cl active | <5 | CORROSIVE TO METALS - Category 1<br>SKIN CORROSION - Category 1B<br>SERIOUS EYE DAMAGE - Category 1<br>SPECIFIC TARGET ORGAN TOXICITY (SINGLE EXPOSURE) (Respiratory tract irritation) - Category 3 |
| sodium hydroxide                       | ≤1 | CORROSIVE TO METALS - Category 1<br>SKIN CORROSION - Category 1A<br>SERIOUS EYE DAMAGE - Category 1                                                                                                 |

### State regulations

**Massachusetts** : The following components are listed: SODIUM HYPOCHLORITE

**New York** : The following components are listed: Sodium hypochlorite

## Section 15. Regulatory information

- New Jersey** : The following components are listed: SODIUM HYPOCHLORITE
- Pennsylvania** : The following components are listed: HYPOCHLOROUS ACID, SODIUM SALT
- California Prop. 65**

This product does not require a Safe Harbor warning under California Prop. 65.

### International regulations

#### Chemical Weapon Convention List Schedules I, II & III Chemicals

Not listed.

#### Montreal Protocol

Not listed.

#### Stockholm Convention on Persistent Organic Pollutants

Not listed.

#### Rotterdam Convention on Prior Informed Consent (PIC)

Not listed.

#### UNECE Aarhus Protocol on POPs and Heavy Metals

Not listed.

### Inventory list

- Australia** : All components are listed or exempted.
- Canada** : All components are listed or exempted.
- China** : All components are listed or exempted.
- Eurasian Economic Union** : **Russian Federation inventory**: All components are listed or exempted.
- Japan** : **Japan inventory (CSCL)**: All components are listed or exempted.  
**Japan inventory (ISHL)**: Not determined.
- New Zealand** : All components are listed or exempted.
- Philippines** : All components are listed or exempted.
- Republic of Korea** : All components are listed or exempted.
- Taiwan** : All components are listed or exempted.
- Thailand** : All components are listed or exempted.
- Turkey** : All components are listed or exempted.
- United States** : All components are active or exempted.
- Viet Nam** : All components are listed or exempted.

## Section 16. Other information

### Hazardous Material Information System (U.S.A.)

|                  |   |   |
|------------------|---|---|
| Health           | / | 0 |
| Flammability     |   | 0 |
| Physical hazards |   | 4 |
|                  |   |   |

Caution: HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. Although HMIS® ratings and the associated label are not required on SDSs or products leaving a facility under 29 CFR 1910.1200, the preparer may choose to provide them. HMIS® ratings are to be used with a fully implemented HMIS® program. HMIS® is a registered trademark and service mark of the American Coatings Association, Inc.

The customer is responsible for determining the PPE code for this material. For more information on HMIS® Personal Protective Equipment (PPE) codes, consult the HMIS® Implementation Manual.

## Section 16. Other information

### National Fire Protection Association (U.S.A.)



### Procedure used to derive the classification

| Classification                                                              | Justification                               |
|-----------------------------------------------------------------------------|---------------------------------------------|
| CORROSIVE TO METALS - Category 1<br>AQUATIC HAZARD (LONG-TERM) - Category 3 | On basis of test data<br>Calculation method |

### History

Date of printing : 30/03/2026

Date of issue/Date of revision : 10/03/2026

Date of previous issue : 17/12/2024

Version : 1

Key to abbreviations : ATE = Acute Toxicity Estimate  
BCF = Bioconcentration Factor  
GHS = Globally Harmonized System of Classification and Labelling of Chemicals  
IATA = International Air Transport Association  
IBC = Intermediate Bulk Container  
IMDG = International Maritime Dangerous Goods  
LogPow = logarithm of the octanol/water partition coefficient  
MARPOL = International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978. ("Marpol" = marine pollution)  
N/A = Not available  
SGG = Segregation Group  
UN = United Nations

References : Not available.

Revision comments : Not available.

Indicates information that has changed from previously issued version.

### Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.