

ABX Pentra CK-MB RTU	A91A00556BEN	20/06/2006
-----------------------------	---------------------	-------------------

1. Product and company identification	
1.1. Identification of the product	
Product name:	ABX Pentra CK-MB RTU
Product code:	Ref.ABX: A11A01643
1.2. Company identification	
	HORIBA ABX - Rue du Caducée - Parc Euromédecine 34184 MONTPELLIER CEDEX 4 - FRANCE Tel: (33) 4 67 14 15 16 Fax: (33) 4 67 14 15 17
1.3. Emergency phone number	
	Contact the nearest first-aid station

2. Composition/information on ingredients	
2.1. Description of the kit	
	<p>REAGENT 1 : Bi-reagent kit :</p> <p>REAGENT 1 : 1 x 20 ml; Imidazole (62.5 mmol/l), Glucose (27.5 mmol/l), N-acetylcysteine (27.5 mmol/l), Magnesium acetate (13.8 mmol/l), EDTA-Na₂ (2.1 mmol/l), NADP (2.75 mmol/l), Hexokinase (≥3kU/l), Polyclonal antibodies (sheep) against human CK-M; inhibiting capacity (≥2.5 kU/l)</p> <p>REAGENT 2 : 1 x 5 ml; Imidazole (160 mmol/l), Creatine phosphate (170 mmol/l), EDTA-Na₂ (2.1 mmol/l), ADP (11 mmol/l), AMP (28 mmol/l), Diadenosine pentaphosphate (56 μmol/l), Glucose-6-phosphate dehydrogenase (≥7.5 kU/l)</p>
2.2. Hazardous ingredients (name, concentration)	
Sodium azide	< 1 g/l for REAGENT 1 and REAGENT 2
2.3. Dangerous preparation according to 67/548/EEC - 99/45/EEC (Yes/No): No	

3. Hazards identification	
Most important hazards:	Dangerous properties cannot be ruled out. However, risks are not to be expected, if the product is handled by qualified and authorized personnel with the necessary precautions for chemicals/diagnostic reagents.
Specific risks	None

4. First-aid-measures	
Inhalation:	No special requirements.
Skin contact:	Wash skin with water.
Eye contact:	Immediately wash eyes, also under the eyelids, with water for at least 15 minutes.
Ingestion:	Give extensive water for drinking. Get medical attention if feeling unwell.

ABX Pentra CK-MB RTU	A91A00556BEN	20/06/2006
-----------------------------	---------------------	-------------------

5. Fire-fighting measures : Product non flammable

Extinguishing media:	Use whatever is required for the surrounding area.
Specific risks:	Not determined.
Special protective equipment:	Not determined.
Additional recommendations:	None.

6. Accidental release measures

Personal precautions:	Wear a lab coat, protective gloves and safety goggles.
Environmental precautions:	Prevent from getting into sewage, water, ground.
Methods for cleaning/absorption:	Dilute with water and absorb with liquid absorbing material. Dispose of according to local regulations. Clean affected area.

7. Handling and storage

7.1. Handling

Technical measures:	No special requirements.
Precautions:	The usual laboratory precautions should be observed.
Handling recommendations:	No special requirements.

7.2. Storage

Technical measures:	No special requirements.
Storage conditions:	Store protected from light and tightly closed at 2-8°C. Do not freeze.
Incompatible materials:	No special requirements.
Packaging materials:	No special requirements.

8. Exposure controls/personal protection

Specific control parameters:	None determined.
Technical measures:	Not determined.
Respiratory protection:	Not required.
Hand protection:	Protective gloves.
Eye protection:	Safety goggles.
Skin protection:	Lab coat.
Hygiene measures:	Change contaminated clothing. Wash hands before and after use.

ABX Pentra CK-MB RTU	A91A00556BEN	20/06/2006
-----------------------------	---------------------	-------------------

9. Physical and chemical properties

Physical state:	Liquid; aqueous clear solution	Autoflammability:	Not determined
Colour:	Colourless	Explosive properties:	Not determined
Odour:	REAGENT 1: Slight smell (Thiol) REAGENT 2: No characteristic odour	Oxidizing properties:	Not determined
pH:	REAGENT 1: 6.0 at 25°C REAGENT 2: 9.0 at 25°C	Vapour pressure:	Not determined
Boiling point:	Not determined	Relative density:	REAGENT 1: 1.0088 g/ml at 20°C REAGENT 2: 1.0388 g/ml at 20°C
Melting point:	Not determined	Solubility:	Not determined
Flashing point:	Not determined	Partition coefficient:	Not determined
Flammability (solid,	Not determined	Other data:	Not determined

10. Stability and reactivity

Conditions to avoid:	Not determined.
Materials to avoid:	Not determined.
Hazardous decomposition products:	Not determined.
Other data:	None.

11. Toxicological information

Acute toxicity:	Not determined.
Chronic toxicity:	Not determined.
Inhalation:	Not determined.
Skin contact:	Not determined.
Eye contact:	Not determined.
Ingestion:	Not determined.
Other data:	Quantitative data about the toxicity of the product are not available. Dangerous properties cannot be ruled out. However, risks are not to be expected, if the product is handled by qualified and authorized personnel with the necessary precautions for chemicals/diagnostic reagents.

12. Ecological information

Mobility:	Not determined.
Degradability:	Not determined.
Bioaccumulative potential:	Not determined.
Ecotoxicology:	Not determined.
Other hazardous effects:	Quantitative data on ecological effects of the product are not available.

ABX Pentra CK-MB RTU	A91A00556BEN	20/06/2006
-----------------------------	---------------------	-------------------

13. Disposal considerations

Surplus or waste (residues):	These reagents contain less than 1g/l of sodium azide as a preservative. As sodium azide may react with lead and copper to form explosive metal azides, this reagent should be disposed of by flushing with copious amounts of water. The product must be disposed of as a laboratory chemical according to local regulations.
Contaminated packaging:	Packaging must be disposed of according to local regulations.

14. Transport information (International regulations)

General information:	None.
By air (IATA):	Not applicable.
By land (ADR):	Not applicable.
By sea (IMDG):	Not applicable.

15. Regulatory information

Symbols:	None	
R phrases:	None	
S phrases:	None	
Substances:		
Other regulatory requirements:	In no way does this information exempt the user from knowing about or applying all the national or international regulations relating to his/her activity.	

16. Other information

This sheet is a complement to instructions for use but does not replace them. Information contained herewith is based on our current knowledge of the product, at the date shown, and is correct to the best of our knowledge.
Furthermore, the user's attention is drawn to the dangers of using this product for anything other than its intended use.
The user must accept the sole responsibility and take precautions accordingly for the use of this product.