

ABX Pentra Low CRP Control	A91A00650CEN	28/01/2008
-----------------------------------	---------------------	-------------------

1. Product and company identification	
1.1. Identification of the product	
Product name:	ABX Pentra Low CRP Control
Product code:	Ref.ABX: A11A01731
1.2. Company identification	
	HORIBA ABX - Rue du Caducée - Parc Euromédecine 34184 MONTPELLIER CEDEX 4 - FRANCE Tel: (33) 4 67 14 15 16 Fax: (33) 4 67 14 15 17
1.3. Emergency phone number	
	Contact the nearest first-aid station.

2. Composition/information on ingredients	
2.1. Description of the kit	
	CONTROL : 4 x 1 ml (ready-to-use): Normal human serum, Human CRP antigen
2.2. Hazardous ingredients (name, concentration)	
	Sodium azide < 0.1%
2.3. Dangerous preparation according to 67/548/EEC - 99/45/EEC (Yes/No): No	

3. Hazards identification	
Most important hazards:	No adverse health effects should occur from routine use of this material in the manner specified by manufacture's instructions. With FDA licensed tests, the human source material in this product has been tested and is not reactive for hepatitis B surface antigen, anti-Hepatitis C, and anti HIV-1/HIV-2. It is recommended to handle this product in accordance with Standard Laboratory Precautions. No hazardous product as specified in Directive 67/548/EEC.
Specific risks:	Skin and/or eye contact, inhalation and ingestion.

4. First-aid-measures	
Inhalation:	Not determined.
Skin contact:	Wash skin with water and soap. In the event of redness or itching, seek medical advice.
Eye contact:	Immediately wash eyes, also under the eyelids, with water for at least 15 minutes. Seek medical advice if irritation persists.
Ingestion:	Wash mouth out with plenty of water and seek medical advice.

ABX Pentra Low CRP Control	A91A00650CEN	28/01/2008
-----------------------------------	---------------------	-------------------

5. Fire-fighting measures

Extinguishing media:	No special requirements.
Specific risks:	Not determined.
Special protective equipment:	Not determined.
Additional recommendations:	None. Thermal decomposition: dangerous decomposition is not anticipated.

6. Accidental release measures

Personal precautions:	No known analysis method can totally guarantee the absence of transmissible pathogenic agent. It is therefore recommended that this product be treated as potentially infectious, and handle observing the usual safety precautions.
Environmental precautions:	Place all spill residue into a suitable container, seal, label and hold for disposal.
Methods for cleaning/absorption:	Absorb and place in the appropriate container for infectious wastes. Wash contaminated area with disinfectant.

7. Handling and storage

7.1. Handling

Technical measures:	No special requirements.
Precautions:	The usual laboratory precautions should be observed.
Handling recommendations:	Wear a lab coat, protective gloves and safety goggles. Avoid splashing or spraying this product's component.

7.2. Storage

Technical measures:	Not determined.
Storage conditions:	Store at 2-10°C. Avoid freezing.
Incompatible materials:	No special requirements.
Packaging materials:	No special requirements.

8. Exposure controls/personal protection

Specific control parameters:	Not determined.
Technical measures:	Follow Standard Laboratory Practices concerning manipulation of human origin products.
Respiratory protection:	Mask.
Hand protection:	Protective gloves.
Eye protection:	Safety goggles.
Skin protection:	Lab coat.
Hygiene measures:	Wash hands before and after use.

ABX Pentra Low CRP Control	A91A00650CEN	28/01/2008
-----------------------------------	---------------------	-------------------

9. Physical and chemical properties

Physical state:	Liquid	Autoflammability:	Not determined
Colour:	Reddish	Explosive properties:	Not determined
Odour:	Odourless	Oxidizing properties:	Not determined
pH:	Not determined	Vapour pressure:	Not determined
Boiling point:	Not determined	Relative density:	Not determined
Melting point:	Not determined	Solubility:	Not determined
Flashing point:	Not determined	Partition coefficient:	Not determined
Flammability (solid, gaz):	Not determined	Other data:	None

10. Stability and reactivity

Conditions to avoid:	Not determined.
Materials to avoid:	Not determined.
Hazardous decomposition products:	None.
Other data:	No hazardous reactions under normal use.

11. Toxicological information

Acute toxicity:	Human serum should be regarded as potentially infectious material.
Chronic toxicity:	Not determined.
Inhalation:	Not determined.
Skin contact:	Not determined.
Eye contact:	Not determined.
Ingestion:	Not determined.
Other data:	There is no evidence of any adverse toxicological problems.

12. Ecological information

Mobility:	Not determined.
Degradability:	Not determined.
Bioaccumulative potential:	Not determined.
Ecotoxicology:	Not determined.
Other hazardous effects:	Not determined.

13. Disposal considerations

Surplus or waste (residues):	These reagents contain less than 0.1% sodium azide as a preservative. As sodium azide may react with lead and copper to form explosive metal azides, these reagents should be disposed of by flushing with copious amounts of water. The products must be disposed of according to local regulations.
Contaminated packaging:	Not determined.

ABX Pentra Low CRP Control	A91A00650CEN	28/01/2008
-----------------------------------	---------------------	-------------------

14. Transport information (International regulations)	
General information:	This product is not subject to current regulations for transportation of hazardous goods.
By air (IATA):	Not applicable.
By land (ADR):	Not applicable.
By sea (IMDG):	Not applicable.

15. Regulatory information	
Symbols:	None
R phrases:	None
S phrases:	None
Substances:	None
Other regulatory requirements:	In no way does this information exempt the user from knowing about or applying all the national or international regulations relating to his/her activity.

16. Other information	
<p>This sheet is a complement to instructions for use but does not replace them. Information contained herewith is based on our current knowledge of the product, at the date shown, and is correct to the best of our knowledge.</p> <p>Furthermore, the user's attention is drawn to the dangers of using this product for anything other than its intended use.</p> <p>The user must accept the sole responsibility and take precautions accordingly for the use of this product.</p>	